


evosuperior has been designed to meet highest requirements for state-of-the-art diagnosis and analysis in left and right heart cath procedures.

The system comprises of recording, measurement and calculation of hemodynamic parameters. It is valued for its brilliance in signal acquisition and signal representation.



## Performance characteristics of evosuperior

### Intuitive operation via extended function keyboard

- One single keyboard for the entire workflow – from measuring to reporting
- Quick access to relevant curve sections by means of button commands
- Editing the oxygen saturation of the relevant heart section by one-button-functionality
- Switching between various measuring intervals, changing the gain of the ECG-signal and sweep speed happens by pushing a single button

### Measurements

- Automatic and simultaneous recording of signal sections
- Manual and automatic measuring
- Manual adjustment of measured sections
- Automatic assignment of pullback curves

In addition to 12 standard measuring positions, further measuring points are configurable. The common measuring algorithms can be freely assigned to the newly generated positions.

### Examination room displays

Next to the real-time signals, the slave monitor displays all measurements jointly with their respective conditions, easy to be monitored by the operator in the examination room.

For reporting purposes the system automatically populates a procedure log during the case, which also allows for manual additions.

### Calculations

Calculated parameters are:

- **Pressure gradients**
- Cardiac-Output
- Vascular Resistance
- Valve Opening Area
- Systemic and Pulmonal Flow
- Shunts
- Bodysurface Area BSA

Algorithms and formulas can be edited by the user; retrospectively changed values will be considered for the calculations.

### ECG Signal Management

The ECG signal is continuously measured in real-time, enabling an immediate reaction to an ischemic condition. Due to the stabilization of the zero position (can be switched on/off), reliable results can be obtained even under rough circumstances.

### Full Disclosure File

All signals are automatically recorded in "full disclosure mode"; therefore the whole data set can be reviewed and analyzed retrospectively after examination.

### Safety concept

The evosuperior is equipped with two separate computer systems. In the unlikely event of a system crash, the secondary system takes over the online monitoring function. This concept guarantees a safe and monitored termination of the examination by displaying ECG and the vital parameters as well as hemodynamic signals.

### Standard Interfaces

HL7 (Health Level 7), optional


- Transfer of patient data from HIS to evosuperior
- Transfer of examination report from evosuperior to HIS

DICOM (Digital Imaging and Communications in Medicine)

- Via DICOM Work List Management (WLM) patient data is transferred to the imaging modality
- Via DICOM Modality Performed Procedure Step (MPPS) the dose information is transferred from the modality to evosuperior and documented in the report
- Via DICOM PACS the report is transferred and stored on the PACS.

# evosuperior

High-end heart cathlab system  
two monitors, two PCs



## Components

- Base Unit (Chassis)
  - Central Interface Unit
  - Power Supply, 230V
  - Isolation Transformer
  - 2 PCs
- 2 x 19" TFT-Monitors
- Patient Headbox
- Function Keyboard DE/EN/FR
- Laserprinter
- Mouse
- UPS; APC-Smart-UPS  
1500W/230V
- Options
  - Vitalmonitor for Cardiac-  
Output, NiBP, SPO<sub>2</sub>
  - (Colour-)Laserprinter
- Cabling
  - Serial connection cable to  
patient headbox
  - 4 carbon fibre cables (radio  
transluscent), length= 90cm
- Standard Accesories
  - Transducer
  - ECG adhesive electrodes
  - 2 pressure cables
  - 1 mounting set for  
2 pressure-transducers
  - ECG patient cable with  
electrode leads

## Technical Specs

<b>Classification acc to European MDD (93/42/EWG):</b>	Class IIb
<b>Patient connection box, Preamplifier</b>	
<b>ECG Amplifier</b>	
Number of amplifiers	9
Sampling rate	500 Hz
<b>IBP Amplifier</b>	
Number of amplifiers	4
<b>Signal Display</b>	
Number of channels	Max. 18
Recording programs	Freely programmable
<b>ECG</b>	
Derivations	I, II, III, aVL, aVR, aVF, V1-V6
Amplification	5, 10, 20 mm/mV
Range of heart rate	20 to 240 beats/min
<b>Invasive blood pressure</b>	
Amplification	10, 25, 50, 100, 200, 400 mmHg
<b>Signal Output</b>	
Analog	8 channels selectable
Trigger	Open collector
<b>Computer OS</b>	Windows 7, 32bit
<b>Monitor</b>	2 x 19" TFT
<b>Printer</b>	Laser monochrom
<b>Applied standards</b>	IEC 60601-1:1988 + A1:1991 + A2:1995 UL 60601-1:2003 IEC 60601-1-1:2000 IEC 60601-1-2:2001+ A1:2004 IEC 60601-2-25:1993 + A1:1999 IEC 60601-2-34:2000 IEC 60601-1-4:1996 + A1:1999 IEC 60601-1-6:2004
<b>Patient Safety</b>	Protection class I, type CF according to IEC 60601-1; Patient connections protected against the effects of defibrillation
<b>Dimensions (h x w x d)</b>	
Basic Unit	710 x 560 x 600 mm
Patient Headbox	265 x 200 x 65 mm
<b>Weight</b>	
Base Unit	ca. 100 kg
Patient Headbox	ca. 0.8 kg
Function Keyboard	ca. 2.0 kg
TFT-Monitor	ca. 10.6 kg
<b>Labeling</b>	CE 0197

Notice: Content may be subject to change